

User Manual
MHS-60
Multi-Sample
Homogenization System
More Intelligent, More Efficient



Cautions

Please be acquainted with the following icons before using this instrument.



Warning Danger (please be careful and note safety)



Key Information (please kindly notice)



Caution Electricity



Mind Your Hands

CONTENTS

目 录



Cautions	01
Contents	02
Announcements	03
I. Introduction	04
II. Unpacking Instructions	04
III. Installation	04
IV. Environment	04
V. Hardware	05
1. Control Terminal	05
2. Transparent Shield	05
3. Homogenizing Probe	05
4. Inner Holder	05
5. Inner Top Cover	06
6. Spider Tube Rack	06
VI. Software	06
1. Main Interface	06
2. Settings Interface	07
3. System Settings	08
VII. Operation	08
1. Manual Mode	08
2. Automatic Mode	08
3. Position Detection	09
4. Date & Time Settings	10
VIII. FAQ & Solutions	10
IX. Service	11
Statement	12

Announcements

Please read this user manual and other attached instructions carefully before using this equipment to ensure proper and safe application. Please keep user manual properly for reference anytime.

This equipment should be positioned in a place which meets laboratory management requirements. Please keep a safe distance from the equipment in running status.

To prevent any accidents and damages to equipment and ensure safe operation,

- please do not use it for other unreasonable purposes;
- please do not disassemble and modify it without authorization;
- please keep it away from inorganic chemistry laboratory;
- please be sure to apply original parts and accessories that has been provided or certified by our company.
- before powering on the equipment, please make sure actual voltage is consistent with the voltage value marked on the nameplate and supplied voltage is stable and grounded.
- please power off immediately after the confirmation of abnormality or malfunction.
- please be careful not to drop the equipment, expose it to physical impact or place heavy objects on it.
- please do not wipe or clean the surface with high-corrosive organic solvents.
- if water, reagents or organic solvents are accidentally spilled on the surface of equipment, please immediately wipe them off with a clean cloth.

To ensure its normal operation, please do not perform a homogenization experiment when there is too little or too much liquid in the homogenizing tube.

Within warranty period, if there is any unartificial malfunctions of the equipment, we promise to immediately provide you with satisfying and convenient services.

If you have any questions, please do not hesitate to call our service hotline +86-020-34086606 or send emails to dt@detelogy.com for further assistance.